

## Subcutaneous immunoglobulin (SCIg) therapy using Springfusor device

1. Load syringe into Springfusor by firmly pushing plunger into the opening of the Springfusor device as shown.



2. To secure the syringe in place, turn the barrel of the syringe at a 90 degrees angle, so that the wings at the base of syringe barrel are sitting securely in the device.



3. Syringe should now be firmly inside the Springfusor as shown and will start delivering the SCIg infusion.



© ASCIA 2017

ASCIA is the peak professional body of clinical immunology and allergy specialists in Australia and New Zealand

### Disclaimer

This document has been adapted with permission from resources developed by Princess Margaret Hospital Immunology Department (Department of Health, Western Australia). It has been peer reviewed by ASCIA members and is based on expert opinion and the available published literature at the time of review. Information contained in this document is not intended to replace medical advice and any questions regarding a medical diagnosis or treatment should be directed to a medical practitioner. Development of this document is not influenced by commercial organisations.