

Subcutaneous immunoglobulin (SCIg) therapy

Preparing the infusion set and inserting HlgH Flo[®] device

The HlgH Flo[®] Device:



a. Preparing the HlgH Flo[®] Device

- a. Attach the tubing from the HlgH Flo[®] device to the tip of the syringe. Check that it is secure.



- b. To remove air bubbles and prime the line, keep the syringe upright and push the plunger of the syringe until the fluid fills the line.



If possible try to prime the line just up to the base of the needle. Keeping the needle "dry" may make needle insertion more comfortable.

b. Inserting the needle

- a. Wipe away the numbing (anaesthetic) cream beforehand (if used)
- b. Wash and dry your hands
- c. Hold the wings of the infusion set and remove the protective cap.
- d. Pinch the skin and insert the needle at a 90 degrees angle.



**PINCHING THE SKIN IS IMPORTANT
TO ENSURE THE NEEDLE GOES
INTO FATTY TISSUE AND
NOT INTO MUSCLE
Remember always pinch an inch!**

- e. Ensure the needle is flush with the skin. Secure the needle with either tape or a clear adhesive dressing.
- f. Tape excess tubing in place, leaving some 'give' in case the tubing is pulled.



© ASCIA 2017

ASCIA is the peak professional body of clinical immunology and allergy specialists in Australia and New Zealand

Disclaimer

This document has been adapted with permission from resources developed by Princess Margaret Hospital Immunology Department (Department of Health, Western Australia). It has been peer reviewed by ASCIA members and is based on expert opinion and the available published literature at the time of review. Information contained in this document is not intended to replace medical advice and any questions regarding a medical diagnosis or treatment should be directed to a medical practitioner. Development of this document is not influenced by commercial organisations.